



Endless Love 

To Guard The New Life Safe



## D60/D70/D80 Patient Monitor



Stock code [300314]

**Product configuration** ("✓" represents standard function, while "●" stands for optional accessory.)

	ECG	SpO <sub>2</sub>	Respiration	Non-invasive blood pressure	Temperature	End-tidal carbon dioxide	Print	Screen size
D60	✓	✓	✓	✓	✓	●	●	8 inches
D70	✓	✓	✓	✓	✓	●	●	10 inches
D80	✓	✓	✓	✓	✓	●	●	12 inches

-  HD screen for clear display
-  Touch (optional) & shuttle dual operation for work efficiency improvement
-  Extra wide viewing angle for observation from different angles
-  Intelligent light control of the screen to automatically adapt to the light environment
-  Fanless design to create a quiet environment for the department
-  Easy to disassemble the plug-in box for easy functional expansion
-  360° alarm lamp, easy to observe alarm
-  Optional mounting brackets, convenient to combine with other medical devices



Large font interface  
Display the measured values of each parameter in an eye-catching bold type, easier for medical staff to identify.

**Accurate and reliable physiological parameter monitoring**

- **SpO<sub>2</sub> monitoring**  
Non-invasive, continuous measurement of pulsating oximetry was used to realize the simultaneous monitoring of blood oxygen parameters in patients.
- **CO<sub>2</sub> monitoring (optional)**  
Provide two measurement methods: mainstream and external bypass.
- **ECG monitoring**  
Professional ECG algorithm, suitable for monitoring the 3-lead (optional configuration) and 5-lead (standard configuration) ECG waveform and heart rate.
- **NIBP monitoring**  
Providing dual overvoltage protection for hardware and software to safe and accurate blood pressure monitoring for children and adults over 3 years of age.

